



PHỤ LỤC SẢN PHẨM

STT	TÊN TRANG THIẾT BỊ Y TẾ	ĐƠN VỊ TÍNH	CHUNG LOẠI SẢN PHẨM	TÊN CƠ SỞ SẢN XUẤT/HÃNG, NƯỚC SẢN XUẤT	TÊN CƠ SỞ CHỦ SỞ HỮU/HÃNG, NƯỚC CHỦ SỞ HỮU	TÊN CƠ SỞ BẢO HÀNH	ĐỊA CHỈ CƠ SỞ BẢO HÀNH
<div style="display: flex; justify-content: space-between;"> <div style="width: 15%;">  </div> <div style="width: 80%; text-align: center;"> <p>DỤNG CỤ PANH, KÌM, KẸP DÙNG TRONG PHẪU THUẬT MẮT</p> </div> <div style="width: 15%; text-align: right;">  </div> </div>							
1	Jaffe Capsule Forcep round grasping tips extremely thin 11 mm long shanks.	Cái	IW-883	Indo Webal Surgicals Pvt.Ltd., Ấn Độ	Indo Webal Surgicals Pvt.Ltd., Ấn Độ	Công ty TNHH Thiết bị và Vật tư y tế Đông Dương	Phòng 606, Tháp A, Tòa nhà Indochina Plaza, 241 Xuân Thủy, phường Dịch Vọng Hậu, quận Cầu Giấy, Hà Nội, Việt Nam
2	Polack Corneal Suturing Forceps very delicate 1mm spread 1x2 teeth 0.12mm.	Cái	IW-884				
3	Harms-colibri suturing forceps very delicate 1x2 teeth, 0.12mm with tying platform.	Cái	IW-885				
4	Gill-Colibri forceps very fine pointed tips with tying platform.	Cái	IW-886				
5	Colibri forceps very delicate 1x2 teeth 0.12mm	Cái	IW-887				
6	Botvin Iris forceps for peripheral iridectomy serrated tips 1x2 teeth	Cái	IW-888				
7	Castroviejo Colibri forceps very delicate 1x2 teeth, 0.12mm	Cái	IW-889				
8	Barraquer colibri-very delicate forceps 0.12mm-1x2 teeth, with typing platform 7mm long	Cái	IW-890				
9	Castroviejo suturing forceps 1x2 teeth with typing platform	Cái	IW-891				
10	Bishop harman tissue forceps 1x2 teeth delicate	Cái	IW-892				
11	Pierse type micro forceps straight skeleton fine Available in: 18-0.5, 19-0.2, 20-0.8mm tip	Cái	IW-893				
12	Bonn forceps 1x2 tip, 0.12mm with typing platform	Cái	IW-894				
13	Hoskin-1 Tissue forceps fine grooved tips, straight & curved	Cái	IW-895				
14	Adson tissue forceps, 1x2 teeth straight, 125mm long delicate	Cái	IW-896				
15	St.Martin corneal suture forceps, straight shaft 1x2 teeth with 6.0mm long typing platform	Cái	IW-897				

30	Green Fixation Forceps multiple teeth , 5.0mm & 10.0 mm wide jaws	Cái	IW-911
31	Sceral Fixation / S.I.C.S. Forceps 1x2 teeth 6 mm platform curved	Cái	IW-912
32	Indo Webal Fixation Forceps 2x3 theeth, straight / angular	Cái	IW-913
33	Wills Hopital Utility Forceps with straight shaft crisscross serrated platform 13.0mm long	Cái	IW-915
34	Bonaccolto Utility Forceps with straight shaft lengthwise serrayed jaws 16.0 mm long 1.7 mm wide with cross serrated tips	Cái	IW-916
35	Moorfeld Conjunctival Suture Forceps straight shaft plain serrated jaws with groove	Cái	IW-917
36	Barraquer Cilia Forceps lightly angled textured jaws provide superior grasping capabilities	Cái	IW-918
37	Beer Cilia Forceps with 4.5 mm long rounded platform	Cái	IW-919
38	Arruga Capsule Forceps 2.5 mm cup curved shaped.	Cái	IW-920
39	Prince Advancement Forceps 10.0mm long jaws with three teeth and sliding lock, right / left	Cái	IW-921
40	Jameson Muscle Forceps with slide lock four teeth child size, right / left	Cái	IW-922
41	Worth Advancement Forceps 11.0mm longitudinal serrated jaws with slliding lock, right / left	Cái	IW-923
42	Berke Ptosis Clamp multy curved shafts with straight blades and longitudinal serrations with sliding lock Available in: 20.0 mm, 27.0 mm	Cái	IW-924
43	Snellen Entropiom Forceps with solid lower plate and open upper plate with locking thumb screw, right / left	Cái	IW-925
44	Ayer Chalazion Forceps with 10mm open upper plate & solid lower with locking thumb screw	Cái	IW-926

11
C
H
H
A
A
10/0

45	Lambert Chalazion Forceps round 8mm I.D	Cái	IW-927
46	Desmarres Chalazion Forceps oval 20 x13mm I.D large	Cái	IW-928
47	Hunt Chalazion Forceps round 12mm I.D	Cái	IW-929
48	Indo Webal Universal Lens Folding forceps to fold the lens perfectly and groove provided in the jaws to hold the lens in the center	Cái	IW-930
49	Shepard Lens Holding Forceps with serrated lower jaws gently curved shaft	Cái	IW-931
50	Shepard Lens Holding Forceps with 3mm lower jaw ,1.5mm upper jaw	Cái	IW-932
51	Blaydes Lens Holding Forceps very delicate narrow jaws angled 45°	Cái	IW-933
52	Dodick Lens Holding Forceps angled smooth lower shaft 12mm & upper shaft 11mm to have a firm grip on iol's	Cái	IW-935
53	Kraff Lens Holding Forceps with 3mm curved tips to grasp to the optic & the haptics.	Cái	IW-936
54	Kratz Lens Holding Forceps very delicate narrow jaws curved	Cái	IW-937
55	Tennant Tying Forceps extra delicate smooth jaws 6mm long tying platform round handle with guide pin	Cái	IW-938
56	Mc pherson Tying Forceps 5.00mm long smooth jaws straight with tying platform, big	Cái	IW-939
57	Mc pherson Tying Forceps 5.00mm long smooth jaws straight with tying platform, small	Cái	IW-939A
58	Mc pherson Tying Forceps curved with smooth jaws, ideal for 8-0 to 11-0 suture with extra delicated 6mm tying platform, big	Cái	IW-940

59	Mc pherson Tying Forceps curved with smooth jaws, ideal for 8 - 0 to 11 - 0 suture with extra delicated 6mm tying flatform, small	Cái	IW-940A
60	Kelmen - Mc pherson Forceps 11.0 mm long platform angeled with smooth jaw ideal for implanting iol & 8 - 0 to 11 - 0 suture, big	Cái	IW-941
61	Kelmen - Mc pherson Forceps 11.0 mm long platform angeled with smooth jaw ideal for implanting iol & 8 - 0 to 11 - 0 suture, small	Cái	IW-941A
62	Jaffe Tying Forceps with smooth jaws ideal for 8-0 to 11 - 0 sutures due to 6.00 mm long platform	Cái	IW-942
63	Mc pherson with smooth jaws 11.00 Mm long platform angled, ideal for implanting IOLs	Cái	IW-943
64	Jaffe Capsulorhexis Forceps very delicate round grasping tips extremely thin 11.00mm long shanks	Cái	IW-944
65	Utrata Capsulorhexis Forceps very delicate triangular grasping tips extremely thin 11.00mm long shanks	Cái	IW-946
66	Utrata Capsulorhexis Forceps very delicate triangular grasping tips extremely thin 11.00mm long shanks round handle with guide pin	Cái	IW-947
67	Nevyas Capsulorhexis Forceps cystotomy shaped tips extremely thin 11.00mm long shanks for penetrating and grasping the anterior capsule	Cái	IW-948
68	Masket Capsulorhexis Forceps very delicate tringular grasping tips extremely thin 11.00 mm long vaulted shanks	Cái	IW-949
69	Inamura Capsulorhexis Forceps very delicate grasping tips and extremely thin 11.00 mm long curveed	Cái	IW-950
70	Perone lasik flap forceps abrasive blasted inner surface.	Cái	IW-953



