

LDHI2 B C

LDHI2

For USA: **CONTENT**
B: N-methylglucamine
 400 mmol/L, pH 9.4;
 lithium lactate 62 mmol/L
C: NAD 62 mmol/L

Lactate Dehydrogenase acc. to IFCC ver.2
cobas c systems
REF 05169330190
 System-ID 05 6607 0

CONTENT 750
 2-8 °C **IVD**
 dialog.roche.com **CE** 0123

LOT 00000000
 2000-01-31

Roche

H317,P261,P272,
 P280,P333+P313,
 P362+P364,P501

Roche Diagnostics GmbH
 Sandhofer Str. 116
 D-68305 Mannheim

Distribution in USA by:
 Roche Diagnostics, Indianapolis, IN
 Made in Germany

06781411001(01)
cobas 04

Pos.1 - 05169330101 - V4




Pos.2 - 05169330102 - V1

LDHI2 **REF** (240)05169330190 For USA: Rx only

LOT (10)00000000
 2000-01-31
 2000-01-01

For USA: Store upright

UDI 

GTIN (01)04015630925056

EU: +800 5505 6606
 USA: +1 800 4282336

Roche 05

Pos.3 - 05169330103 - V5

Nhãn sản phẩm

LDHI2 (05169330190)

**Thuốc thử xét nghiệm định lượng lactate
dehydrogenase**

LDHI2

Mã sản phẩm: 05169330190

Số lưu hành:

Chủ sở hữu số lưu hành và nhập khẩu:


Công ty TNHH Roche Việt Nam


Tầng 27, Tòa nhà Pearl Plaza, 561A Điện Biên Phủ,
Phường 25, Quận Bình Thạnh, Tp. Hồ Chí Minh

Chủ sở hữu: Roche Diagnostics GmbH, Sandhofer
Strasse 116, 68305 Mannheim, Đức

Nhà sản xuất: Roche Diagnostics GmbH, Sandhofer
Strasse 116, 68305 Mannheim, Đức

Xuất xứ: Đức

Hướng dẫn bảo quản: xem  , **Số lô:** xem ,

Hạn sử dụng: xem  (năm-tháng-ngày / năm-tháng)
trên nhãn

Hướng dẫn sử dụng xem tại

Các thông tin khác xem trong hướng dẫn sử dụng

Hỗ trợ kỹ thuật vui lòng liên hệ 1800599991