

HOW TO USE BREATHING SYSTEM FILTERS

- Prevent viruses, bacteria, and cross-infection from patient to another when using ventilator, and respiratory function meter.
- Medical plastic outer material; The inner membrane prevents viruses and bacteria 99.999%.
- For use with ventilators: One end is attached to the air line (breathing line), one end is attached to the endotracheal tube or connector.
- Used for respiratory function meter: One end is attached to the body of the device, and the other end is for the patient to suck on to blow.

GIÁM ĐỐC



Nguyễn Thị Huệ Chi