

MẪU NHÃN


Tên cơ sở đăng ký lưu hành: **Công ty TNHH Thương mại - Dịch vụ Kỹ thuật Lục Tỉnh**
Địa chỉ: 849 Trần Xuân Soạn, Phường Tân Hưng, Quận 7, TP. HCM

Ngày 18 tháng 05 năm 2022

NỘI DUNG NHÃN BẰNG TIẾNG VIỆT

IVD kiểm tra pha loãng
QuikLYTE Integrated Multisensor Dilution Check

Số lưu hành: xxxxx

Số lô: Xem  trên nhãn gốc

HSD: Xem  trên nhãn gốc

Bảo quản: Xem  trên nhãn gốc

Chủ sở hữu: Xem  trên nhãn gốc

Nước chủ sở hữu: Mỹ

Nước xuất xứ: Mỹ

Chủ sở hữu số lưu hành: Công ty TNHH Thương mại - Dịch vụ Kỹ thuật Lục Tỉnh
849 Trần Xuân Soạn, Phường Tân Hưng, Quận 7, TP. HCM

SĐT: 028 37755740

Cảnh báo, Hướng dẫn sử dụng, Hướng dẫn bảo quản xin xem trên Hướng dẫn sử dụng tại:

www.xxxxx.xxx

Người đại diện hợp pháp của cơ sở
Phó giám đốc




Tống Thị Bích Tuyên

CTN_L_QuikLYTE Dilution Check_DM
Part Number: 10866592_05 (755640.101)
Drawing Number: L Carton
Template Number: NA
Manufacturing Site: GLA

COLORS:
Black
PMS 321
PMS3258

10866592_05


L-01

SIEMENS

QuikLYTE®
Dilution Check



SIEMENS


Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA
siemens-healthineers.com
EC REP
Siemens Healthcare Diagnostics
Manufacturing Ltd.
Chapel Lane
Swords, Co. Dublin, Ireland

REF S640
SMN 10444872

SIEMENS

QuikLYTE® Integrated Multisensor
Dilution Check



siemens.com/eifu

IVD

2°C-8°C



CONTENTS

1 x 50 mL

For US **RxOnly**
Intended Use: Used to monitor and adjust for the dilution ratio on Dimension® clinical chemistry systems with the QuikLYTE® module.
Reagents:
Na+ 140 mmol/L
K+ 4 mmol/L
Dimension®, QuikLYTE® (USA)

Rev. G

SIEMENS
QuikLYTE® Integrated Multisensor
Dilution Check

SIEMENS

REF **S640**

LOT



QuikLYTE®
Dilution Check



[siemens.com/eifu](https://www.siemens.com/eifu)

IVD



2°C — 8°C

CONTENTS

50 mL

10866591_06




Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA
[siemens-healthineers.com](https://www.siemens-healthineers.com)

LABEL_BTL_4.25_QuikLYTE Dilution Check_DM
Part Number: 10866591_06 (755640.001)
Dimensions: 4.25" x 2.5"
Template Number: NA
Manufacturing Site: GLA

COLORS:
Black
PMS 321
PMS 3258

Product Name	Dimension QuikLYTE Dilution Check
SMN	10444872
REF	S640
Label	UDI Label (Template: 11275239)



Disclosure and Confidential Statement

This document contains information which is strictly confidential and proprietary. It is to be used exclusively for registration by Siemens Regulatory Representatives, Local Agency Reviewers and Health Authorities. The distribution, copying or publication of any information in this document is prohibited.