

MẪU NHÃN

Tên cơ sở đăng ký lưu hành: **Công ty TNHH Thương mại - Dịch vụ Kỹ thuật Lạc Tĩnh**
Địa chỉ: 849 Trần Xuân Soạn, Phường Tân Hưng, Quận 7, TP. HCM

Ngày 16 tháng 05 năm 2022

NỘI DUNG NHÃN BẰNG TIẾNG VIỆT

**IVD kiểm tra hệ thống
Dimension Flex reagent cartridge CHK**

Số lưu hành: xxxxx

Số lô: Xem  trên nhãn gốc

HSD: Xem  trên nhãn gốc

Bảo quản: Xem  trên nhãn gốc

Chủ sở hữu: Xem  trên nhãn gốc

Nước chủ sở hữu: Mỹ

Nước xuất xứ: Mỹ

Chủ sở hữu số lưu hành: **Công ty TNHH Thương mại - Dịch vụ Kỹ thuật Lạc Tĩnh**
849 Trần Xuân Soạn, Phường Tân Hưng, Quận 7, TP. HCM

SĐT: 028 37755740

Cảnh báo, Hướng dẫn sử dụng, Hướng dẫn bảo quản xin xem trên Hướng dẫn sử dụng tại:
www.xxxxx.xxx

Người đại diện hợp pháp của cơ sở
Phó giám đốc



Tổng Thị Bích Tuyên

CHK

REF **DF179**
SMN 10481507

SIEMENS

SIEMENS

Dimension® clinical chemistry system
Flex® reagent cartridge

 siemens.com/eifu

IVD

CONTENTS 8



WARNING
H317
P280,
P302+P352,
P333+P313,
P362+P364



+H505CHK2X

CHK
2°C 8°C

For US **RxOnly**
CHK Flex® reagent cartridge
Intended Use: For checking proper functioning of the fluid metering and photometric measuring system on the Dimension® clinical chemistry system. It is run daily as part of the System Check routine.
Reagents: Red dye with buffer and preservative.

Dimension® and Flex® (USA)

SIEMENS

 Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 USA
siemens-healthineers.com

EC REP
Siemens Healthcare Diagnostics
Manufacturing Ltd.
Chapel Lane
Swords, Co. Dublin, Ireland

CTN_8UP_FLEX_CHK_DM
Part Number: 10864891_02 (717006.179)
Drawing Number: D Carton (10916373)
Template Number: NA
Manufacturing Site: GLA

COLORS:
Black
PMS 185
PMS 321

Rev. H

SIEMENS

Dimension® clinical chemistry system
Flex® reagent cartridge

CHK



10864891_02

CHK Flex[®]

IVD

LOT

XXXXXX

2°C/8°C

Σ 20

Siemens Healthcare Diagnostics Inc.
Newark, DE 19714 USA



XXXX-XX-XX

12345678



Dimension® (USA), Flex® (USA)

Dimension[®] clinical chemistry system
Flex[®] reagent cartridge

IVD

CHK

2°C/8°C

LOT

XXXXXX



XXXX-XX-XX




WARNING

Siemens Healthcare Diagnostics Inc.
Newark, DE 19714 USA

11554551_01

Product Name	Dimension CHK Flex reagent cartridge
SMN	10481507
REF	DF179
Label	UDI Label (Template: 10916782)



Disclosure and Confidential Statement

This document contains information which is strictly confidential and proprietary. It is to be used exclusively for registration by Siemens Regulatory Representatives, Local Agency Reviewers and Health Authorities. The distribution, copying or publication of any information in this document is prohibited.