

ALT

B
ALT

C

For USA: **CONTENT**
B: TRIS 224 mmol/L, pH 7.3;
L-alanine 1120 mmol/L;
albumin (bovine) 0.25 %
LDH (microb.) ≥ 45 µkat/L
C: 2-Oxoglutarate 94 mmol/L;
NADH ≥ 1.7 mmol/L

Alanine Aminotransferase
cobas c systems

REF 05850797190
System-ID 05 7501 0

CONTENT 1100

2-8 °C **IVD**

 dialog.roche.com

COBAS C and COBAS
are trademarks of Roche.



Roche Diagnostics GmbH
Sandhofer Str. 116
D-68305 Mannheim

Distribution in USA by:
Roche Diagnostics, Indianapolis, IN
Made in Germany

cobas[®]

0443202001(3)

LOT 00000000
2000-01-31



ALT

REF (240)05850797190

For USA: Rx only

LOT (10)00000000

2000-01-31

2000-01-01

UDI



GTIN (01)04015630928927



For USA: Store upright

EU: +800 5505 6606
USA: +1 800 4282336



05

Nhãn sản phẩm

ALT (05850797190)

**Thuốc thử xét nghiệm định lượng alanine
aminotransferase**

ALT

Mã sản phẩm: 05850797190

Số lưu hành:

Chủ sở hữu số lưu hành và nhập khẩu:


Công ty TNHH Roche Việt Nam


Tầng 27, Tòa nhà Pearl Plaza, 561A Điện Biên Phủ,
Phường 25, Quận Bình Thạnh, Tp. Hồ Chí Minh

Chủ sở hữu: Roche Diagnostics GmbH, Sandhofer
Strasse 116, 68305 Mannheim, Đức

Nhà sản xuất: Roche Diagnostics GmbH, Sandhofer
Strasse 116, 68305 Mannheim, Đức

Xuất xứ: Đức

Hướng dẫn bảo quản: xem  , **Số lô:** xem **LOT** ,

Hạn sử dụng: xem  (năm-tháng-ngày / năm-tháng)
trên nhãn

Hướng dẫn sử dụng xem tại

Các thông tin khác xem trong hướng dẫn sử dụng

Hỗ trợ kỹ thuật vui lòng liên hệ 1800599991