



Document Information:

Document Author: Maschat, Cordula {DOMC~Mannheim}
Business Area / Unit: Roche Professional Diagnostics
Confidentiality: Confidential
Document Class: Device Master Record
Document Type: Device Master Record
Document Creator: Maschat, Cordula {DOMC~Mannheim}
Document Lifecycle Status: Signed
Valid From: 01-Jun-2021 05:47:27 (UTC)
Valid To:
Document Title: C_0201_RSP_03507343000
Document Number: 0000000000001008018000535
Document Version: 03
Template: No

Global Group: Regulatory Affairs
Global SubGroup: Specification Prod. (Registr.)
Local Group: Reg. Specification
Site: RDG Germany
Department: C_DOMCC
Document Description: Tina-quant IgA Gen.2 (IGA-2) - "Only for use in Regulatory Affairs"

Electronic Signatures:

Signed By: maschac1 (Cordula Maschat {DOMCSD})
Role: Reviewer
Signature Differentiation: Quality Control
Signed Date: 01-Jun-2021 04:33:14 (UTC)

Signed By: kneisses (Sandra Kneissel {DOMC})
Role: Approver
Signature Differentiation: Quality Control
Signed Date: 01-Jun-2021 05:39:33 (UTC)



Specification

Cat. No.: 03507343190

Product: **Tina-quant IgA Gen.2 (IGA-2)**

Parameter	Specification
Reagents-working solutions	
Buffer solution (B):	
pH-value (25 °C):	7.90 - 8.10
Anti-human IgA antibody/buffer solution (C):	
pH-value (25 °C):	7.90 - 8.10

Parameter**Specification**

Performance test**Recovery on cobas c analyzer(s) Standard application:**

Precinorm Protein (Cat. No. 10557897122): 90 - 110 %

Precipath Protein (Cat. No. 11333127122): 90 - 110 %

Recovery on cobas c analyzer(s) Sensitive application:

Precinorm Protein (Cat. No. 10557897122): 90 - 110 %

Precipath PUC (Cat. No. 03121291122): 90 - 110 %

Also required:

Calibrator f.a.s. Proteins (Cat. No. 11355279216)

**Roche Diagnostics GmbH
Quality Control Mannheim**Sandhofer Strasse 116
D-68305 Mannheim
Phone +49-621-759-0