
















## FITTINGS FOR COLD LIGHT CABLES

<b>Fittings Endoscope</b>	<b>Connection</b>	<b>Ordering no.</b>
<i>RUDOLF / Storz / Olympus / Winter &amp; Ibe / Aesculap</i>		<b>IM220-020</b>
<i>Wolf / Henke-Sass / Dyonics / ID = 9.0 mm</i>		<b>IM220-021</b>
<i>ACMI / British Standard / Codman / ID = 7.7 mm</i>		<b>IM220-022</b>
<i>Schott / Zeiss</i>		<b>IM220-024</b>
<i>Stryker</i>		<b>IM220-025</b>

## FITTINGS FOR LIGHT SOURCES

<b>Fittings Light Source</b>	<b>Connection</b>	<b>Ordering no.</b>
<i>RUDOLF / Storz / Aesculap</i>		<b>IM220-040</b>
<i>Wolf / Henke-Sass / Dyonics / Comeg / Medicon</i>		<b>IM220-041</b>
<i>ACMI / British Standard / Codman / Welch Allyn / Stryker</i>		<b>IM220-042</b>
<i>Olympus</i>		<b>IM220-043</b>
<i>Winter &amp; Ibe</i>		<b>IM220-046</b>
<i>Pentax</i>		<b>IM220-047</b>
<i>Fuji</i>		<b>IM220-048</b>
<i>Volpi</i>		<b>IM220-049</b>