

## Specification

### Handheld Base Unit (08408793190)

Specification	Handheld Base Unit
Height	82 mm / 3.22 in (max.)
Width	120 mm / 4.72 in (max.)
Length	228 mm / 8.97 in (max.) 295 mm / 11.61 in (max.) with extension
Weight	340 g / 11.99 oz (max.) 410 g / 14.46 oz (max.) with extension
User interface	LED
Interfaces	Charging contacts IR port RJ45 Ethernet USB type B
Operating temperature	12 to 32 °C / 54 to 90 °F
Storage temperature (long-term)	-5 to +45 °C / 23 to 113 °F
Humidity (operating)	10 to 85% (no condensation)
Relative humidity (storage)	< 95%
Input voltage	12 V DC
Input current	1.25 A (max.)
Data transfer rate	IR: 9.6-115.2 Kbit/s Ethernet: 10/100 Mbit/s (auto-negotiate) full-duplex USB: 480 Mbit/s