

# Elecsys Anti-Tg

**cobas e analyzers**

REF 07026919190      CONTENT  $\Sigma$  300

COBAS, COBAS E and ELECSYS are trademarks of Roche.

R1 19.7 mL  
R2 19.7 mL  
M 14.1 mL

Roche Diagnostics GmbH  
Sandhofer Str. 116  
D-68305 Mannheim

Distribution in USA by:  
Roche Diagnostics,  
Indianapolis, IN

Made in Germany      001

IVD CE 2-8 °C

<http://e-labdoc.roche.com>

Pos.1

# ATG

REF 07026919190

GTIN 04015630939619

LOT 00000000

2000-01-31

2000-01-01

03

Pos.2

# Elecsys Anti-Tg

REF 07026919190

For USA: CONTENT

R1 Biotinylated Tg (human) 0.200 mg/L  
R2 Monoclonal anti-Tg antibodies (human) labeled with ruthenium complex 0.620 mg/L  
M Streptavidin-coated microparticles 0.72 mg/mL

Rx only      001

Pos.3

# ATG

07026919190

R1  
R2  
M

IVD

LOT 00000000

2000-01-31

01

Pos.4

# ATG

1310065

$\Sigma$  300

2-8 °C

01

Pos.5

## Nhãn sản phẩm

**Elecsys Anti-Tg (07026919190)**

### Thuốc thử xét nghiệm định lượng kháng thể kháng thyroglobulin

#### Elecsys Anti-Tg

**Mã sản phẩm:** 07026919190

**Số lưu hành:**

**Chủ sở hữu số lưu hành và nhập khẩu:**


Công ty TNHH Roche Việt Nam


Tầng 27, Tòa nhà Pearl Plaza, 561A Điện Biên Phủ,  
Phường 25, Quận Bình Thạnh, Tp. Hồ Chí Minh

**Chủ sở hữu:** Roche Diagnostics GmbH, Sandhofer  
Strasse 116, 68305 Mannheim, Đức

**Nhà sản xuất:** Roche Diagnostics GmbH, Sandhofer  
Strasse 116, 68305 Mannheim, Đức

**Xuất xứ:** Đức

**Hướng dẫn bảo quản:** xem  , **Số lô:** xem  ,

**Hạn sử dụng:** xem  (năm-tháng-ngày / năm-tháng)  
trên nhãn

Hướng dẫn sử dụng xem tại

Các thông tin khác xem trong hướng dẫn sử dụng

Hỗ trợ kỹ thuật vui lòng liên hệ 1800599991