



a member of the Roche group

QUALITY CONTROL SPECIFICATION

Product name: VENTANA HE 600 Hematoxylin

Ventana REF: 07024282001

Roche Number: 07024282001

Functional Staining on
small intestine

The chromatin (hetero and euchromatin) must be visible and well-defined (crisp).

The nuclear membrane must be visible and contiguous.

The nuclei bluish-purple to purple in color.

The intensity must be strong enough to allow clear demonstration of nuclear detail at 10X power, but not so strong as to cause loss of chromatin granularity or excessive background staining.


Mucin staining must be present and must be lighter than the nuclei.

At least 3 of the 4 small intestine sections must meet these criteria.

pH

2.45 – 2.54


Vincent A. Faddis
Quality Department Manager


Date