



Certificate of Analysis

Abbott Ireland Diagnostics Division,
Lisnamuck,
Longford,
Co. Longford,
Ireland

Product Name: Cholesterol2
List number: 04S9230
Lot number: 36397UD00
Reagent Lot Number: 36397UD00
Expiry Date: 2023-05-10

SCP / revision 4S920030/004

Component name	List Number	Lot Number	Production date	Expiry date
Chol2 R1	04S92R	36396UD00	2022-02-10	2023-05-10

Acceptance Table: Grand Mean Concentration (mg/dL)

Description	Lower Limit	Upper Limit	Value	Result
REF Control 1	235.30	254.70	240.34	Acceptable
REF Control 2	102.70	110.98	104.45	Acceptable
REF_GC1_Lv5	580.08	641.14	604.14	Acceptable

This product has been manufactured, including packaging and quality control in full compliance with ISO 13485, IVD Directive 98/79/EC and GMP guidelines as applicable.

Note: The specifications in this certificate reflect the internal Abbott release specifications. When performing the assay in a laboratory the specifications listed in the package insert apply.

Quality:

Kathleen Greer
Printed name and signature

Date:

01 Mar 22