



Technical Specification

Product name: Access WASH BUFFER II

Part number: A16792

Intended use: Access Wash Buffer II is intended for use with the Access Immunoassay Systems and specific Access immunoassay reagents.

Storage condition: 15-30°C

Parameter	Specification
General Inspection	Label information must be correct
Appearance	Liquid, Colorless to pale yellow, Clear, Odorless
pH	8.11 - 8.69
Bioburden	≤ 100 CFU/mL
Volume	$\geq 1,950$ mL/bottle
Package	4 bottles/box

Ekta Ghimire

Analyst II Regulatory Affairs

Beckman Coulter, Inc.