





Reagent Kit II





Electrolyte Analyzer Supporting Reagent (Ion Selective Electrode Method and Pressure Method)

Intended For URIT-910CPlus Electrolyte Analyzer

Reagent Composition:		  
URIT A Standard Solution A	400mL/Pack	
URIT B Standard Solution B	100mL/Pack	
URIT F Electrode Cleaning Solution	20mL/Pack	



Wear disposable gloves to avoid contact with potentially infectious material.

URIT Medical Electronic Co.,Ltd.

http://www.urit.com

Rev:06/2018-C1



MẪU NHÃN PHỤ

Tên trang thiết bị y tế: Thuốc thử định lượng điện giải Natri, Kali, Chlorua, Calci, pH

Chủng loại:

Reagent Kit II: Standard Solution A: 400 ml/pack

Standard Solution B: 100 ml/pack

Electrode cleaning solution: 20 mL/pack

Số đăng ký lưu hành:

Tên và địa chỉ chủ sở hữu số lưu hành, cơ sở bảo hành:

CÔNG TY TNHH ĐẦU TƯ VÀ THƯƠNG MẠI NGUYỄN HƯNG

Số 37, nhà B, Tổ 25 Tập thể VH Nghệ Thuật, Phường Mai Dịch, Quận Cầu Giấy, Thành phố Hà Nội.

Nhà sản xuất: **URIT Medical Electronic Co., Ltd**

Địa chỉ: No.D-07 Information Industry District, High-Tech Zone, Guilin, Guangxi 541004, Trung Quốc

Xuất xứ: Trung Quốc

Số lô: Xem **LOT**

Ngày sản xuất: Xem

Hạn sử dụng:

Thông tin cảnh báo, hướng dẫn sử dụng, hướng dẫn bảo quản: Xin xem hướng dẫn sử dụng kèm theo