
Quality Control Specification

Product name: CEA Mouse Monoclonal Antibody

Ventana PN : 760-4594

Roche Number: 6433316001

*Date completed: 04/02/12 supercedes N/A

Test	Acceptance Criteria
Functional Staining performed on: Colon	Intensity difference between sample and a reference must be less than or equal to 0.5.
	Background difference between sample and a reference must be less than or equal to 0.5.
	Cytoplasmic staining pattern observed.
Visual Inspection	Transparent liquid and free of foreign material.
Bioburden	Visible bacterial colonies must be less than 100CFU/ml.


Quality Assurance Representative

02 APRIL 2012
Date