



TEST REPORT

HORIBA, Ltd. certifies the items mentioned below.

Customer: FURUNO ELECTRIC Co., Ltd

1. Specification

Product Name : ISE Cleaner 3200295053 3200504560
 Lot No. : 220111
 Production date : 2022/1/11
 Expiration date : 2023/1/11

2. Inspection Results

pH and Osmotic pressure measurement Good
 pH Within pH 8.0 ±0.2
 Osmotic Pressure Within 330 ± 20mOsm.

3. Overall Judgment Passed

4. Other matters

Room temp: 21.6 °C Relative Humidity: 35 %

Date: 2022-01-18

Inspected by: Taishi Karadka

Approved by: Akemi Mori

HORIBA, Ltd.