

---

## Quality Control Specification

---


**Product name:** Uroplakin III Rabbit Monoclonal Antibody

**Ventana PN :** 760-4533

**Roche Number:** 6419232001

\*Date completed: 10/09/12 supercedes \_\_\_\_\_

Test	Acceptance Criteria
Functional Staining performed on: Transitional Cell Carcinoma	Intensity difference between sample and a reference must be less than or equal to 0.5.
	Background difference between sample and a reference must be less than or equal to 0.5.
	Membranous, Cytoplasmic staining pattern observed.
Visual Inspection	Transparent liquid and free of foreign material.
Bioburden	Visible bacterial colonies must be less than 100CFU/ml.

  
\_\_\_\_\_  
Quality Assurance Representative

Oct. 09, 2012  
Date