

Tên thiết bị:

Ống thông động mạch

Số lưu hành:**Chủ sở hữu số lưu hành:**

Công ty TNHH TM Sixmurs HTH Việt Nam,
Số nhà 9, ngõ 7, đường Lê Đức Thọ,
Phường Mỹ Đình 2, Quận Nam Từ Liêm,
TP. Hà Nội, Việt Nam

Cơ sở sản xuất:

Arrow Internacional de Chihuahua S.A. de
C.V.,
Ave. Washington 3701, Interior Circuito
Industrial Alta Tecnologia Edificios
2, 4, 36 And 40, Colonia Panamericana
CHIHUAHUA, Chihuahua Mexico
CP 31200

Chủ sở hữu sản phẩm:

Arrow International LLC (subsidiary of
Teleflex Incorporated),
3015 Carrington Mill Blvd., Morrisville, North
Carolina 27560, Mỹ


Xuất xứ:

Mexico

Số lô:

Xem **LOT** trên nhãn gốc

Hạn sử dụng:

Xem  trên nhãn gốc

Cơ sở bảo hành:

Không áp dụng

**Thông tin cảnh báo, hướng dẫn sử dụng,
hướng dẫn bảo quản:**

Xem tờ Hướng dẫn sử dụng kèm theo sản
phẩm



GIÁM ĐỐC

Phạm Thị Thu Hằng

- CS Katetrizační set pro přístup přes arteria radialis
- DA Kateteriseringssæt til arteria radialis
- NL Arteria radialis katheterisatieset
- FR Jeu de cathétérisme pour l'artère radiale
- DE Besteck zur Katheterisierung der Arteria radialis
- EL Σετ καθετηριασμού κερκδικής αρτηρίας
- HU Katéterező felszerelés arteria radialis-hoz
- IT Set per cateterismo dell'arteria radiale
- PL Zestaw do cewnikowania tętnicy promieniowej
- PT Conjunto de Cateterização da Artéria Radial
- RU Комплект для катетеризации лучевой артерии
- SK Sada pre radiálnu arteriálnu katetrizáciu
- SL Komplet za arterijsko radialno katetrizacijo
- ES Conjunto de cateterización de arteria radial
- SV Set för radiell artärkatetrisering
- TR Radyal Arter Kateterizasyon Seti

- CS Katetrizační set pro přístup přes arteria radialis
- DA Kateteriseringssæt til arteria radialis
- NL Arteria radialis katheterisatieset
- FR Jeu de cathétérisme pour l'artère radiale
- DE Besteck zur Katheterisierung der Arteria radialis
- EL Σετ καθετηριασμού κερκδικής αρτηρίας
- HU Katéterező felszerelés arteria radialis-hoz
- IT Set per cateterismo dell'arteria radiale
- PL Zestaw do cewnikowania tętnicy promieniowej
- PT Conjunto de Cateterização da Artéria Radial
- RU Комплект для катетеризации лучевой артерии
- SK Sada pre radiálnu arteriálnu katetrizáciu
- SL Komplet za arterijsko radialno katetrizacijo
- ES Conjunto de cateterización de arteria radial
- SV Set för radiell artärkatetrisering
- TR Radyal Arter Kateterizasyon Seti

REF RA-04020 Radial Artery Catheterization Set

ARROW®

QTY: 25

Contents:

Catheters: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque Polyurethane over
22 Ga. TW Introducer Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide
Wing Clips: Suture

Rx only.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Contraindications: Inadequate collateral arterial circulation.

Sterile. Single use — do not resterilize. Sterilized by ethylene oxide. Does not contain natural rubber latex. Fluid path components are non-pyrogenic. Contains no medication.

 Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

 Teleflex Medical
IDA Business and
Technology Park
Athlone, Ireland



IA04020-1380 (2021-04)

REF RA-04020 Radial Artery Catheterization Set

ARROW®

QTY: 25

Contents:

Catheters: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque Polyurethane over
22 Ga. TW Introducer Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide
Wing Clips: Suture

Rx only.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Contraindications: Inadequate collateral arterial circulation.

Sterile. Single use — do not resterilize. Sterilized by ethylene oxide. Does not contain natural rubber latex. Fluid path components are non-pyrogenic. Contains no medication.

 Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

 Teleflex Medical
IDA Business and
Technology Park
Athlone, Ireland



IA04020-1380 (2021-04)