

DO NOT BEND**ARROW**REF **SAC-01222** (IPNXXXXXX)**22** Ga.**12** cm catheter length**Arterial Catheterization Set****Contents:**

- 1: Indwelling Catheter: 22 Ga. x 12 cm Radiopaque Polyurethane
 1: Spring-Wire Guide, Marked: .021" (0.53 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on Both Ends)
 1: Introducer Needle: 22 Ga. x 1-1/2" (4 cm) XTW

Catheter Priming Volume: 0.18 mL Priming Volume is approximate.

Not made with natural rubber latex.
 Fluid path components are non-pyrogenic.

Contraindications: Inadequate collateral arterial circulation.

- cs** Arteriální katetrizační set
da Arterielt kateteriseringsæt
nl Arteriële katheterisatieset
et Arterikanüüli komplekt
fi Valtimokatetrisarja
fr Jeu de cathétérisme artériel
de Besteck zur Katheterisierung von Arterien
el Σετ αρτηριακού καθετηριασμού
hu Artériás katéter készlet
it Set per cateterismo arterioso
lt Arterijų kateterizavimo rinkinys
no Arteriell kateteriseringssett
pl Zestaw do cewnikowania tętnic
pt Conjunto de cateterização arterial
ro Set de cateterizare arterială
ru Набор для артериальной катетеризации
sk Arteriálna katetrizačná sada
sl Pribor za arterijsko kateterizacijo
es Conjunto de cateterización arterial
sv Artärkateteriseringsset
tr Arteriyel Kateterizasyon Seti



(01)XXXXXXXXXXXXXXXXXX
 (17)YYMMDD
 (10)LotNumber

LOT LotNumber

YYYY-MM-DD

Arrow International LLC
 Subsidiary of Teleflex Incorporated
 3015 Carrington Mill Blvd.
 Morrisville, NC 27560 USA

Teleflex[®]
 www.teleflex.com
 LBL062165 (2021-06)

EC REP Teleflex Medical
 IDA Business and
 Technology Park
 Dublin Road, Athlone
 Co. Westmeath, Ireland
 2797

REF **SAC-01222****Arterial Catheterization Set****22** Ga.**12** cm catheter length**DO NOT BEND****ARROW**