



Document Information:

Document Author: Schnabel, Helena {DOMC~Mannheim}
Business Area / Unit: Roche Professional Diagnostics
Confidentiality: Confidential
Document Class: Device Master Record
Document Type: Device Master Record
Document Creator: Schnabel, Helena {DOMC~Mannheim}
Document Lifecycle Status: Signed
Valid From: 02-Jun-2021 11:32:05 (UTC)
Valid To:
Document Title: C_0201_RSP_06749992000
Document Number: 0000000000001200000027206
Document Version: 02
Template: No

Global Group: Regulatory Affairs
Global SubGroup: Specification Prod. (Registr.)
Local Group: Reg. Specification
Site: RDG Germany
Department: C_DOMCC
Document Description: Tina-quant Lambda Gen.2 LAMB2 - "Only for use in Regulatory Affairs"

Electronic Signatures:

Signed By: schnabeh (Helena Schnabel {DOMCSD})
Role: Reviewer
Signature Differentiation: Quality Control
Signed Date: 02-Jun-2021 11:23:04 (UTC)

Signed By: kneisses (Sandra Kneissel {DOMC})
Role: Approver
Signature Differentiation: Quality Control
Signed Date: 02-Jun-2021 11:25:26 (UTC)



Specification

Cat. No.: 06749992190

Product: **Tina-quant Lambda Gen.2 (LAMB2)**

Parameter	Specification
Reagents-working solutions	
Buffer (B):	
pH-value (25 °C):	7.8 - 8.2
Antibodies in buffer (C):	
pH-value (25 °C):	7.3 - 7.7

Parameter	Specification
-----------	---------------

Performance test

Recovery on COBAS INTEGRA analyzer(s):

Precinorm Protein (Cat. No.: 10557897122): 90 - 110 %

Precipath Protein (Cat. No.: 11333127122): 90 - 110 %

Recovery on cobas c analyzer(s):

Precinorm Protein (Cat. No.: 10557897122): 90 - 110 %

Precipath Protein (Cat. No.: 11333127122): 90 - 110 %

Also required:

C.f.a.s. Proteins (Cat. No. 11355279216)

**Roche Diagnostics GmbH
Quality Control Mannheim**

Sandhofer Strasse 116
D-68305 Mannheim
Phone +49-621-759-0