



Certificate of Analysis

Abbott Ireland Diagnostics Division,
Lisnamuck,
Longford,
Co. Longford,
Ireland

Product Name: Prealbumin
List number: 01E0224
Lot number: 39417UD00 **SCP / revision** 1E020024/003
Reagent Lot Number: 39417UD00
Expiry Date: 2023-01-17

Component name	List Number	Lot Number	Production date	Expiry date
Palb R1	01E02R	39418UD00	2021-12-17	2023-01-17
Palb R2	01E02T	39419UD00	2021-12-17	2023-01-17

Acceptance Table: Grand Mean Concentration (mg/dL)

Description	Lower Limit	Upper Limit	Value	Result
REF Control 1	17.16	20.58	18.89	Acceptable
REF Control 2	22.08	26.48	24.33	Acceptable
REF Control 3	30.77	36.91	34.09	Acceptable
REF Control 4*	47.25	57.75	52.71	Acceptable

*Grand Mean REF_Control_4 Concentration acceptance criteria is $\pm 10.00\%$ from the Control 4 Expected Value as per SCP 1E020024


This product has been manufactured, including packaging and quality control in full compliance with ISO 13485, IVD Directive 98/79/EC and GMP guidelines as applicable.

Note: The specifications in this certificate reflect the internal Abbott release specifications. When performing the assay in a laboratory the specifications listed in the package insert apply.

Quality:


 Printed name and signature

Date:


 15 Apr 22