

Certificate of Analysis

Product Name: TRANSFERRIN CLINICAL CHEMISTRY	Abbott List Number: 01E0422
OEM Part Number: 01E0422	Abbott Commodity Number: 01E0422
OEM Lot Number: 29327UN22	Expiration Date: 2023-04-11
Storage: 2 - 8°C	Package Insert No./Rev. G95903 / R03
Qty Manufactured: 46 Kits	

Specification	Specification Limit	Result
Barcode Scan	Acceptable verification on c8000	Complies Pass
Bio-Rad L1	138 - 168 mg/dL	155 Pass
Bio-Rad L1	Lot Number	68971
Bio-Rad L2	248 - 303 mg/dL	275 Pass
Bio-Rad L2	Lot Number	68972
Bio-Rad L3	335 - 410 mg/dL	368 Pass
Bio-Rad L3	Lot Number	68973
Control 4: % Difference	-10 - 10 %	0 Pass
Micro: Beg Sample (R1)	<30 CFU/mL	0 Pass
Micro: End Sample (R1)	<30 CFU/mL	0 Pass
Micro: Beg Sample (R2)	<30 CFU/mL	0 Pass
Micro: End Sample (R2)	<30 CFU/mL	0 Pass

For in-vitro diagnostic use.

This lot meets all current release specifications per quality document MS.MASTER342.

The above product has been manufactured in compliance with Quality Systems Regulations, passed all Thermo Fisher final release specifications, and has met all Abbott acceptance requirements.

This product contains animal source material. The caprine material used in the manufacture of this product is derived from goat antiserum collected from U.S. sourced goats.

This product contains R1 and R2 components. The component name, lot number, and expiration date matches the kit name, lot number, and expiration date listed on this COA.



Quality Representative

Date: 2022-05-06