

ISE Cleaner

Intended Use

ISE Cleaner is ISE Wash solution on BX series clinical chemistry analyzers.

Instruction

Please follow by Operation Manual for cleaning procedure.
ISE Cleaner should be handled by trained personnel.

Composition

NaCl < 1.0%
Protease < 0.5%

Storage & Shelf Life after First Opening

Store in original container protected from direct sunlight in a dry, cool dark place and well-ventilated area, away from incompatible materials and food and drink. Keep container tightly closed and sealed until ready for use. Containers that have been opened must be carefully resealed and kept upright to prevent leakage. Do not store unlabeled containers. Use appropriate containment to avoid environmental contamination.

Store at temperature between 5°C to 30°C.

After opening, this product is stable for 6 months.

Shelf life: 12 months after manufacturing date

Warnings & Precautions

Before use, please read the operator's manual of BX series analyzers carefully.

Put on appropriate personal protective equipment. Do not ingest.

No eating, drinking, and smoking in working areas. Workers should wash hands and face before eating, drinking, and smoking. Remove contaminated clothing and protective equipment before entering eating areas.

Additional Special Equipment

ISE Cleaner is intended for use with the BX-series clinical chemistry analyzers.

Manufacturer

SYSMEX CORPORATION



1-5-1, Wakino-hama-Kaigandori, Chuo-ku, Kobe 651-0073, Japan

Manufacturing Site:

HORIBA, Ltd.

2 Miyano-higashi, Kisshoin, Miamo-ku, Kyoto 601-8510, Japan

Authorised Representative/Distributors

Singapore: **SYSMEX ASIA PACIFIC PTE. LTD**
9 Tampines Grande, #06-18, Singapore 528735

Japan: **SYSMEX CORPORATION**
1-5-1 Wakino-hama-Kaigandori, Chuo-Ku, Kobe 651-0073, Japan

Product Information

CB322736 ISE Cleaner 100 mL × 1

Issued date or revision

10/2021

This PI's information is based on the following sources:

- Safety Data Sheet: 07/2015
- Product label