

## KNIFE



*Other Names:* Perforator, Separator

*Category:* Knife

*Used:* Knife instrument is a specialized device that surgeons use to pry up the muscles in order to expose the bone and perform osteotomy.

*Description:* Knife instrument is a specialized device that surgeons use to pry up the muscles in order to expose the bone and perform osteotomy. It includes the technical characteristics: Sharp Tip, Curved, Blade Working Surface, Hollow Grip Handle, Stainless Steel, Reusable.

*Notes:* Inspect instruments before each use for broken, cracked, chipped or worn parts. Never allow blood to dry on the instruments. Rinse the instruments with distilled, filtered water. Never use tap water as it contains minerals that could leave a residue (stain) on the instrument surface.