

ANTRUM CURETTE



Other Names: currett

Category: Curettes

Used: Curette instrument used to scrape away tissue or bone.

Description: Curette instrument is a specialized device that orthopedic surgeons use to gouge out abnormal bone growths and other hard tissues, in order to perform certain surgical procedures. Circular & Oval tips to suit multiple approaches, 2.0 to 8.0mm spoon sizes to scrape different surfaces.

Notes: Inspect instruments before each use for broken, cracked, chipped or worn parts. Never allow blood to dry on the instruments. Rinse the instruments with distilled, filtered water. Never use tap water as it contains minerals that could leave a residue (stain) on the instrument surface.