



JEOL Ltd., : 3-1-2 Musashino, Akishima Tokyo 196-8558 JAPAN
Tel: (+81)42-542-2303 Fax: (+81)42-542-3132

Lamp Coolant C

Intended Use

Lamp coolant containing anti rust agent.

Instruction of Use

Pour the lamp coolant C to the lamp coolant tank of the instrument according to the Instruction manual.

Product Information

Part number : 780656661

Pack size : 400mL

Store at 5 to 25°C

Composition

Propylene glycol < 10%

Anti-rust solution < 0.01%

Preservative < 0.05%

Storage

Store in a dark, cool place and avoid direct sunlight.

Keep the lid tightly closed.

Store at 5 to 25°C

Disposal

Dispose of residual product with plenty of water.

If substantial volume needs to be disposed, dispose it through approve waste disposal companies.

Manufacturer

JEOL Ltd.

3-1-2 Musashino, Akishima, Tokyo 196-8558, Japan

Revision date

19 / Dec / 2017