
Quality Control Specification

Product name: MUM1 (MRQ-43) Rabbit Monoclonal Antibody

Ventana PN: 760-4529

Roche Number: 06419208001

*Date completed: 8/4/2017 supercedes 4/2/2012

Test	Acceptance Criteria
Functional Staining performed on: Diffuse Large B-cell Lymphoma, Stomach	Intensity difference between sample and a reference must be less than or equal to 0.5.
	Background difference between sample and a reference must be less than or equal to 0.5.
	Nuclear, cytoplasmic staining pattern observed.
Visual Inspection	Transparent liquid and free of foreign material.
Bioburden	Visible bacterial colonies must be less than 100CFU/ml.
ELISA	0.5 µg/ml – 5.0 µg/ml



Quality Assurance Representative

8-29-17
Date

