

**SANYRENE**  
**(Corpitolinol 60)**

**PREVENTION OF PRESSURE ULCERS AND IMPAIRED SKIN INTEGRITY**

**□ DESCRIPTION:**

Oil composed of Corpitolinol 60: hyperoxygenated glycerides of essential fatty acids (linoleic acid (60%), tocopherol (vitamin E)): 99%.

Anise perfume: 1%.

Dosing pump bottle.

**□ PROPERTIES:**

**SANYRENE** restores the protective hydrolipidic film of the skin and accelerates the renewal of epidermal cells.

**SANYRENE** strengthens, soothes and protects the areas of skin which are at risk.

**SANYRENE** rubs in gently without pain, it has a prolonged action.

Hypoallergenic.

**□ INDICATIONS:**

Prevention of pressure sore.

Treatment and prevention of blanching erythema (reversible redness) of the skin, caused by pressure, friction and shear.

Treatment and prevention of dry and dehydrated skins, skin fragility due to maceration conditions, skins with low cutaneous resistance and trophic impairment of skin (i.e. Dry and scaly dermatosis, surrounding wound skin of vascular ulcers).

**□ INSTRUCTIONS FOR USE:**

The use of **SANYRENE** does not prevent from using other therapies or strategies for the prevention of pressure sores.

**SANYRENE** should be used as soon as there is a risk of pressure sores developing, of blanching erythema or of any of the indications above.

Apply one spray of **SANYRENE** (Corpitolinol 60) on average on the area of skin which is at risk. Rub the substance in gently for one minute using the finger tips. Repeat this procedure on each area at risk 3 to 4 times a day, at each position change or as often as needed.

**□ WARNINGS AND PRECAUTION FOR USE:**

Do not use in case of skin damage.

Do not use if the bottle is broken.

External use. Do not swallow.

**□ CONTRA-INDICATION:**

Known sensitivity to the product.



**TECHNICAL DATA:**

NATURE OF CONTROLS	SPECIFICATIONS
Characters	Amber colored, slightly viscous oil, with a characteristic anise odor
Peroxide value	50 to 150 mEq/Kg
Acid value	≤ 1.0
Shelf life	3 years

**REGULATORY STATUS:**

CE Marked Medical device.

**SANYRENE** has been designed and is manufactured by:  
URGO Laboratories – 21300 CHENOVE - FRANCE