


MASTER ALARM PANEL MODEL MPJ



OVERVIEW


MPJ type Master alarm panel is a fully integrated multiplexed system which is capable of 24 hours monitoring for all of medical gas services. In the event of abnormal supply conditions, its alarm lamp illuminates and initiates a buzzer warning.

■ NAMES AND FUNCTIONS



No.	Name	No.	Name
1	Box	5	Alarm confirmation button
2	Face plate	6	Test button
3	Gas name (1~8 gases)	7	Power lamp
4	Alarm lamp	8	Normal lamp

■ SEQUENCE FLOW CHART



* Specification and appearance shall be subject to change without notice.



Head office : 2-3-16 Nishikanda, Chiyoda-ku, Tokyo, Japan 101-0065 /TEL : +81-3-3556-8122 /FAX : +81-3-3556-1572

International department : 5-1-1Maigaoka, Kokuraminami-ku, Kitakyushu, Fukuoka, Japan 802-0823 /TEL+81-93-963-8331/FAX : +81-93-963-8291

Branches : Hokkaido/Tohoku/Tokyo/Nagoya/Osaka/Hiroshima/Kyushu

Factory : Kokura(Fukuoka)

URL : <http://www.central-uni.co.jp>