

Medical Gas Monitoring System

AREA MONITOR APJ



OVERVIEW


Medical gas area monitor is a fully integrated multiplexed system which is capable of 24 hours monitoring for all of medical gas services. In the event of abnormal supply conditions, its alarm lamp illuminates and initiates a buzzer warning.

NOTE


- 1~8 gases are available as shown on below illustrations S type for 1~4gas and L type for 5~8gas.
- The display indicates the pressure condition in the pipeline HIGH or LOW
- Each gas alarm segment is distinguished by a different color on the display. Green for oxygen, blue for nitrous oxide, yellow for medical air, black for vacuum, gray for nitrogen, brown for surgical tool air, etc.. Other color distinctions for specific reasons are possible.

NAMES AND FUNCTIONS

S type (1~4gas)



L type (5~8gas)




No.	Name	No.	Name
1	Power lamp	5	Gas name
2	Normal lamp	6	Pressure display window (digital)
3	Silence button	7	Abnormal lamp (high/low, overload)
4	Test button		

DISPLAY

Pressure value



Abnormal lamp



* Specification and appearance shall be subject to change without notice.

C&U CENTRAL UNI CO.,LTD.

Head office : 2-3-16 Nishikanda, Chiyoda-ku, Tokyo, Japan 101-0065 /TEL : +81-3-3556-8122 /FAX : +81-3-3556-1572

Oversea division : 5-1-1Maigaoka, Kokuraminami-ku, Kitakyushu, Fukuoka, Japan 802-0823 /TEL+81-93-963-8331/FAX : +81-93-963-8291

Branches : Hokkaido/Tohoku/Tokyo/Nagoya/Osaka/Hiroshima/Kyushu

Factory : Kokura(Fukuoka)

URL : <http://www.central-uni.co.jp>

C&U