

INSTRUCTION OF USE

- Make sure all instruments is sterilized before using in hospital .
- Using Instruments in according to the product's recommended purpose, complying with technical specifications to avoid warping and malfunction.
- After soaking in chemicals, the instruments must be rinsed with clean water and dried.
- Instruments can be autoclave or doing low temperature disinfection With drying stage.
- Attention: For long-term use of tools, do not soak chemicals for too long.
- Storage: Keep Instruments in a dry place, away from high moisture of wet places.

Thanks for choosing ADI's Products!

ADI INDUSTRY
Partner

HEAD OFFICE.

Shahab Pura Road Opp. Gepeco Coloni
51310, Sialkot Pakistan
Ph: +92-52-325 93 83
Fax: +92-52-325 93 83
Cell:+92-321-6140483
Email: Info@adiindustry.com
URL: www.adiindustry.com

PRODUCTION UNIT.

Industrial Estate 51310, Sialkot, Pakistan.
Ph: +92-52-325 93 83
Fax: +92-52-325 93 83
Cell:+92-321-6140483
Email: Info@adiindustry.com
URL: www.adiindustry.com

BANKERS.

Habib Bank Limited
Address: Industrial Estate Branch.zip. 51310
Sialkot, Pakistan. Swift Code (kjhatyre)
JS Bank
Address: Industrial Estate Branch.zip. 51310
Sialkot, Pakistan. Swift Code (kjhatyre)