

STERIS



**PROLYSTICA® 2X
CONCENTRATE NEUTRAL DETERGENT
DIRECTIONS FOR USE**

DIRECTIONS FOR USE

Prolystica 2X Concentrate Neutral Detergent may be manually or automatically dispensed at 1/8-1/2 oz/gal (1-4 mL/L) depending on water quality and application.

Medical devices should be sterilized or disinfected by appropriate means in a separate subsequent step.

STORAGE AND DISPOSAL

Store in a cool place, below 32°C (90° F). Keep out of the reach of children. Dispose in accordance with local, state and federal regulations. Empty containers should be rinsed and discarded in a waste container or offered for recycling if possible. Do not use empty containers for other purposes.

SERVICE

Sales

Service is one of the most important ways to verify consistent quality of the facility's performance and operation. A tailored service program by STERIS provides effective, trouble-free operations.

Technical

Contact STERIS for information regarding product applications, the product in general or copies of the Material Safety Data Sheet (MSDS). For information regarding human and environmental exposure, call STERIS's St. Louis operations regulatory affairs group at 314-535-1395.

For further information, please contact:

STERIS®



STERIS Corporation
5960 Heisley Road
Mentor, OH 44060-1834 • USA
440-354-2600 • 800-548-4873
www.steris.com