

# SRI TRANG GLOVES (THAILAND) PUBLIC COMPANY LIMITED

10 Soi 10 Phetkasem Road, Hat Yai, Hat Yai, Songkhla 90110 Thailand  
Tel. (66) 074-344-663 Fax. (66) 074-344-677

SCT.PD.SP.05.065-140222 R.03  
PRODUCT SPECIFICATION SHEET : NCCOGF-E-OS-G-NS

DAR-SCT-2022-0299  
PAGE 1 of 1

## GLOVE SPECIFICATION

**Product Specification Group Code :** NCCOGF-E-OS-G-NS  
**Product Code :** DNCFEOG (With Color)  
**Product name :** Nitrile Gloves, Powder Free, Polymer Coated, Non Sterile  
**Grade :** OS Spec, General Purpose Grade

### GENERAL INFORMATION

**Type :** Single-use disposable or protective glove  
**Marking :** Gloves are not marked to designated size  
**Shape :** Ambidextrous  
**Material :** Nitrile Butadiene Rubber  
**Cuff :** Rolled cuff  
**Surface :** Textured at Finger Tip

<b>Thickness :</b> (single wall according to ASTM D6319, ASTM D3767)	<b>Finger :</b> Min. 0.05mm	<b>Palm :</b> Min. 0.05mm
---	-----------------------------	---------------------------

<b>Length :</b> (according to ASTM D6319, ASTM D3767)	XS, S	: Min. 220mm
	M, L, XL	: Min. 230mm

<b>Tensile Strength :</b> (according to ASTM D6319, ASTM D412, ASTM D573)	Before Aging	: Min. 14MPa
	After Aging	: Min. 14MPa

<b>% Elongation :</b> (according to ASTM D6319, ASTM D412, ASTM D573)	Before Aging	: Min. 500%
	After Aging	: Min. 400%

<b>Available Sizes &amp; Palm width (mm) :</b> XS : 70+/-10, S : 80+/-10, M : 95 +/-10, L : 110+/-10, XL : 120+/-10 (according to ASTM D6319, ASTM D3767)
--

**Storage:** 3 years storage life (date of expiration see on the box). Keep area cool, dry and dust free, avoid ventilation and storage close to photocopy equipment. Copper-ions discolor the glove. Protect gloves against ultraviolet light sources, as sunlight and oxidizing agents. Storage above 30°C will lead to accelerated aging and should be avoided.

### STANDARDS, GUIDELINES AND QUALITY CERTIFICATION

(Every mentioned standard is used in the latest edition.)

**Certification :** ISO 9001, ISO 13485

**Conformity to standards :** ASTM D6319 (Dimension and Physical properties)

Please report any incidents with the products to [STGT\\_TA@srtranggroup.com](mailto:STGT_TA@srtranggroup.com) and [regulatory@srtranggroup.com](mailto:regulatory@srtranggroup.com) without undue delay. Our Product Safety Officers will review all information provided and might contact the original informant for further clarification.

# SRI TRANG GLOVES (THAILAND) PUBLIC COMPANY LIMITED

10 Soi 10 Phetkasem Road, Hat Yai, Hat Yai, Songkhla 90110 Thailand  
Tel. (66) 074-344-663 Fax. (66) 074-344-677

SCT.PD.SP.05.065-140222 R.03

PRODUCT SPECIFICATION SHEET : NCCOGF-E-OS-G-NS

DAR-SCT-2022-0299

Appendix of specification

## APPENDIX : 2

**Product Specification Code :** NCCOGF-E-BLK-OS-GG-NS  
**Product Code :** DNCFEOG BLACK  
**Material :** Nitrile Butadiene Rubber  
**Treatment and Feature :** Powder Free, Polymer Coated, Offline Chlorinated

## SAMPLING, INSPECTION AND FINAL RELEASE INFORMATION

Product Specification Code	Glove color	Characteristics			Test-Method
		Barrier Defects*	Major Visual Defects*	Minor Defects*	
NCCOGF-E-BLK-OS-GG-NS	BLACK	AQL ≤ 10	AQL ≤ 10	-	ASTM D5151

### Barrier Defects :

Inspection level G I (ISO 2859-1)

Holes that results in the malfunction of glove as a protection barrier either physically or not physically affect the appearance.

### Major Visual Defects :

Inspection level G I (ISO 2859-1)

Defined defect that disrupt the glove's function or defect that cause an uncomfortable feeling to the user.

### Minor Defects :

Inspection level G I (ISO 2859-1)

Visual defects aggregated are non-conformities of lower degree of concern, which do not prevent correct use of gloves. Minor cosmetic defect which does not affect on the glove's performance

Remark :

\* Defect Characteristics following In-house practice.