





Product Specification

Name of goods	Interactive ADL recognition rehabilitation												
Model name	MOTOCOQ	Manufacturer, Country	CyberMedic Corp. / Korea										
<p><input type="checkbox"/> Use</p> <p>The system is capable of ADL training, cognitive training and hand function training through the exercise load using tools commonly used in everyday life</p> <p><input type="checkbox"/> Features</p> <ul style="list-style-type: none"> - Interactive training programs such as real life using tangible training tools. - Ten(10) kind of detachable tools(circle(L, S size), cylindrical, ball shape, letter of L, air tube, steering wheel, pinch, button, power grip, and touch screen). - ADL Training function through motion recognition by tool and exercise resistance control. - Three kind of training(ten hand function course, Thirty cognitive training course, six ADL course) - Set program difficulty(1~5 level) - Show the evaluation result of tools and item by item.(hand-eye coordination, active resistance, attention, information processing, classification, memory, recognition, sense of direction, language) - User data management - Upper assistance using five axis arm support (option) - Height control using electrical table(option) 	   												
<p><input type="checkbox"/> Product specifications</p> <ol style="list-style-type: none"> 1) Dimension : 350(L)x340(W)x400(H)mm 2) Power of adapter : DC24V, 1.5A 3) Communication type : more than USB 2.0 3) Ten(10) kind of detachable tools(circle(L, S size), cylindrical, ball shape, letter of L, air tube, steering wheel, pinch, button, power grip, and touch screen) 4) Three kind of training(ten hand function course, Thirty cognitive training course, six ADL course) 	<p><input type="checkbox"/> Product configuration</p> <table border="0"> <tr> <td>1) Main body</td> <td style="text-align: right;">1ea</td> </tr> <tr> <td>2) Ten tools</td> <td style="text-align: right;">1set</td> </tr> <tr> <td>3) 21 inches all-in-one touch PC</td> <td style="text-align: right;">1ea</td> </tr> <tr> <td>4) 24V adapter</td> <td style="text-align: right;">1ea</td> </tr> <tr> <td>5) USB cable</td> <td style="text-align: right;">1ea</td> </tr> <tr> <td>6) User manual</td> <td style="text-align: right;">1ea</td> </tr> </table> <p><input type="checkbox"/> The others</p> <ol style="list-style-type: none"> 1) Warranty : 1 years after installation 2) Installation and education : Free <p style="text-align: right;">(Only domestic)</p>	1) Main body	1ea	2) Ten tools	1set	3) 21 inches all-in-one touch PC	1ea	4) 24V adapter	1ea	5) USB cable	1ea	6) User manual	1ea
1) Main body	1ea												
2) Ten tools	1set												
3) 21 inches all-in-one touch PC	1ea												
4) 24V adapter	1ea												
5) USB cable	1ea												
6) User manual	1ea												