



Document Information:

Document Author: Limburger, Marion {DQXP~Penzberg}
Business Area / Unit: Roche Professional Diagnostics
Confidentiality: Confidential
Document Class: Device Master Record
Document Type: Device Master Record
Document Creator: Limburger, Marion {DQXP~Penzberg}
Document Lifecycle Status: Signed
Valid From: 15-Oct-2020 06:54:17 (UTC)
Valid To:
Document Title: C_08_RSP_07027915190_EN
Document Number: 0000000000001200000156593
Document Version: 03
Template: No

Global Group: Regulatory Affairs
Global SubGroup: Specification Prod. (Registr.)
Local Group: Reg. Specification
Language: English
Site: RDG Germany
Department: C_DOMS
Document Applies To: C_Standardization/ QC Penzberg
Document Description: Elecsys Testosterone II - "Only for use in Regulatory Affairs"

Electronic Signatures:

Signed By: limburgm (Marion Limburger {DQXPFB})
Role: Author
Signature Differentiation: Quality Assurance
Signed Date: 14-Oct-2020 15:06:16 (UTC)

Signed By: wickerts (Sabine Wickert {DOMSW})
Role: Reviewer
Signed Date: 14-Oct-2020 17:17:52 (UTC)

Signed By: tassonir (Richard Tassoni {DOMS})
Role: Approver
Signed Date: 15-Oct-2020 06:46:54 (UTC)

SPECIFICATION
for
Elecsys Testosterone II

0702 7907 190
0702 7915 190
0894 6370 190
0520 0067 190
0894 6353 190

Part 1 (M)

Appearance	sediment with clear supernatant
Streptavidin-coated microparticles	648 - 792 µg/ml

Part 2 (R1)

Appearance	clear solution
Biotinylated monoclonal anti-testosterone antibody	proved by performance test
pH-value	5.8 – 6.2

Part 3 (R2)

Appearance	clear solution
Testosterone derivative ~ruthenium complex	proved by performance test
pH-value	5.8 – 6.2

Performance test

Recovery of Testosterone in
(based on assigned value)

PreciControl Universal Level 1	87 – 113 %
PreciControl Universal Level 2	87 – 113 %