



## Document Information:

Document Author: Goetz, Martina {DOMS~Penzberg}  
Business Area / Unit: Roche Professional Diagnostics  
Confidentiality: Confidential  
Document Class: Device Master Record  
Document Type: Device Master Record  
Document Creator: Ioannidis, Ioannis {DQXP~Penzberg}  
Document Lifecycle Status: Signed  
Valid From: 12-Feb-2020 13:17:55 (UTC)  
Valid To:  
Document Title: C\_08\_RSP\_07027737190\_EN  
Document Number: 0000000000001200000173296  
Document Version: 04  
Template: No

Global Group: Regulatory Affairs  
Global SubGroup: Specification Prod. (Registr.)  
Local Group: Reg. Specification  
Language: English  
Site: RDG Germany  
Department: C\_DOMS  
Document Applies To: C\_Standardization/ QC Penzberg  
Document Description: Elecsys Prolactin II - "Only for use in Regulatory Affairs"

## Electronic Signatures:

Signed By: wickerts (Sabine Wickert {DOMSW})  
Role: Author  
Signed Date: 12-Feb-2020 12:43:02 (UTC)

Signed By: ioannidi (Ioannis Ioannidis {DQXPFB})  
Role: Reviewer  
Signature Differentiation: Quality Assurance  
Signed Date: 12-Feb-2020 13:05:56 (UTC)

Signed By: tassonir (Richard Tassoni {DOMS})  
Role: Approver  
Signed Date: 12-Feb-2020 13:09:59 (UTC)

# SPECIFICATION

for

**Elecsys Prolactin II**

---

0702 7737 190 / 0702 7737 214

0702 7729 190

**0320 3093 190** / 0320 3093 214

## Part 1 (M)

Appearance	sediment with clear supernatant
Streptavidin-coated microparticles	648 - 792 µg/mL

## Part 2 (R1)

Appearance	clear solution
Biotinylated monoclonal anti-prolactin antibody	proved by performance test
pH value	6.5 - 7.5

## Part 3 (R2)

Appearance	clear solution
Monoclonal anti-prolactin antibody labeled with ruthenium complex	proved by performance test
pH value	6.5 - 7.5

## Performance test

Recovery of Prolactin in  
(based on assigned value)

Control 1	
PreciControl Universal Level 1	84 - 116 %
PreciControl Universal Level 2	84 - 116 %