



Document Information:

Document Author:	Schwaiger, Nicole {DOMS~Penzberg}
Business Area / Unit:	Roche Professional Diagnostics
Confidentiality:	Confidential
Document Class:	Device Master Record
Document Type:	Device Master Record
Document Creator:	Schwaiger, Nicole {DOMS~Penzberg}
Document Lifecycle Status:	Signed
Valid From:	30-May-2022 06:58:55 (UTC)
Valid To:	
Document Title:	C_08_RSP_03277356190_EN
Document Number:	0000000000001008018000271
Document Version:	04
Template:	No
Global Group:	Regulatory Affairs
Global SubGroup:	Specification Prod. (Registr.)
Local Group:	Reg. Specification
Language:	English
Site:	RDG Germany
Department:	C_DOMS
Document Applies To:	C_Standardization/ QC Penzberg
Document Description:	Prolactin II CalSet - "Only for use in Regulatory Affairs"

Electronic Signatures:

Signed By:	schwaig (Nicole Schwaiger {DOMSTO})
Role:	Author
Signature Differentiation:	Quality Control
Signed Date:	30-May-2022 04:02:47 (UTC)
Signed By:	wickerts (Sabine Wickert {DOMSW})
Role:	Reviewer
Signed Date:	30-May-2022 05:45:51 (UTC)
Signed By:	tassonir (Richard Tassoni {DOMS})
Role:	Approver
Signed Date:	30-May-2022 06:47:20 (UTC)

SPECIFICATION

for

Prolactin II CalSet

0327 7356 190

Part 1 (PRL Cal 1)

Appearance

yophilized substance

Prolactin

approx. 2 μ IU/mL

Part 2 (PRL Cal 2)

Appearance

yophilized substance

Prolactin

approx. 2000 μ IU/mL

cobas e 801 analyzer: The exact lot-specific calibrator values are encoded in the electronic barcode and available via the cobas link.

All other analyzers: The exact lot-specific calibrator values are encoded in the barcode as well as printed on the enclosed (or electronically available) calibrator barcode sheet.