

Quality Control

Certificate of Analysis

QC release date: 2021-09-23

Product: Cholesterol FS
Product Code: 11300

Cat. No.:

Delivery of kits: No of kits:
Kit Lot/Los:
Kit Expiry:

Delivery of bulk: Volume:
Processing Lot:

Results for:

Reagent Lot No.: 31330
Expiry: 2023-09-24

Examination:

Appearance: Reagent : clear, pink
Odour: Reagent : slight smell (Phenol)
pH: Reagent (25 °C ± 1 °C): 6.71


Absorbance: Reagent : 0.010
546 nm/ 15-25°C

Time to reach endpoint ≤ 10 minutes at 37 °C

Control sera: TruLab N TruLab P
Lot: 29695 Lot: 29696

Fresh (at 37°C) 125 mg/dL 186 mg/dL

Target values 126 108 - 143 188 161 - 214



Additional remarks: _____

This lot meets our specifications.
This certificate of analysis is a computer printout and has therefore not been signed
by hand.

Sign.: Dr. G. Gorka