

This document contains confidential, proprietary information of Ventana Medical Systems, Inc. It may not be copied or reproduced without prior written permission from Ventana Medical Systems, Inc.	<b>Document Part No.</b>	<b>RFT-006-01</b>	<b>REV</b>	<b>Page</b>
	<b>RA Parent Doc No.</b>	<b>RI-006</b>	<b>A</b>	<b>1 of 1</b>

**TITLE: D134717 Rev 2 – Quality Control Specification**

---

**Product Name** Reaction Buffer Concentrate (10X)  
**Ventana**  REF 950-300; 951-300  
**Roche Number** 05353955001; 05353955214

Final Acceptance Testing pH acceptable range must be within 7.6 ± 0.1 at 22°C  
 Conductivity acceptable range must be within 19.07 mS/cm ± 2.4

Visual Inspection The reagent must be clear and free of turbidity.  
  
 The reagent must be free of foreign material and substantially free of particulate matter.

---

**Quality Department Manager**

**Date**

< End of Document >