

HOW TO USE WIRE GUIDED STONE EXTRACTION BALLOON

- Used in endoscopy to remove stones from the bile duct.
- Gravel traction ball, with 2 contrast bands at the beginning and end of the ball, with 3 channels (guidewire channel, balloon pump channel, and contrast solution injection channel).
- Ball pump diameter: 15-18-20mm
- Location of contrast agent: above the balloon
- Compatible with guidewire goes inside the tool
- There is a marking at the end of the tool (V-marking), connection to the scope (C-hook), designed in the shape of a V (V-Sheath).

GIÁM ĐỐC



The stamp is circular with a red border. Inside the border, the text "CÔNG TY TNHH MỘT THÀNH VIÊN HUE CHI GROUP" is written around the top and sides. In the center, there is a red logo of a hand holding a flame. Below the logo, the text "HUE CHI GROUP" is written. At the bottom of the stamp, the text "MSDN: 0307712254" is visible. A blue ink signature is written over the stamp.

Nguyễn Thị Huệ Chi