

HOW TO USE DISPOSABLE SUCTION CONNECTING TUBE

- Used to connect to a suction machine to aspirate fluid during surgery.
- The product has been sterilized, do not use it when the packaging is opened. Use only one time.
- Tear off the packaging.
- Connect to a suction machine and suction head.
- Use as directed by the doctor.

GIÁM ĐỐC



The stamp is circular with a red border. Inside the border, the text "CÔNG TY TNHH MỘT THÀNH VIÊN" is written along the top arc, and "HUE CHIGROUP" is written along the bottom arc. In the center of the stamp is a red logo featuring a stylized flame or bird shape. Below the logo, the text "MSDN: 0307712254" is visible. A blue ink signature is written over the stamp, extending to the right.

Nguyễn Thị Huệ Chi