



## Closed Wound Drainage System (Spring Type)


### Introduction

**GIÁM ĐỐC**  
*Phạm Công Toàn*

This product includes 3-spring evacuator, PVC tubing, Y connector, PVC drainage tube and stainless steel trocar.

The design of 3-spring evacuator minimizes the Collapse-Resistance.

The PVC drainage tube inside can be replaced by our other silicone drains, if needed.



Product	Size	Sterile	Ref. code
Closed wound drainage system	200 ml + 10FR	√	HTJ0101S
	200 ml + 12FR	√	HTJ0102S
	200 ml + 14FR	√	HTJ0103S
	200 ml + 16FR	√	HTJ0104S
	200 ml + 18FR	√	HTJ0105S
	400 ml + 10FR	√	HTJ0201S
	400 ml + 12FR	√	HTJ0202S
	400 ml + 14FR	√	HTJ0203S
	400 ml + 16FR	√	HTJ0204S
	400 ml + 18FR	√	HTJ0205S