

FORMAT

This document contains confidential, proprietary information of Ventana Medical Systems, Inc. It may not be copied or reproduced without prior written permission from Ventana Medical Systems, Inc.

Document Part No.	RFT-001-02	REV	Page
Parent Doc No.	RI-001	B	1 of 1

TITLE: D127963 Rev. B - QUALITY CONTROL SPECIFICATION

Product name: Protease 2

Ventana REF : 760-2019

Roche Number: 05266696001

Absorbance

Reference and Sample Mean must score in the range of 0.100 and 0.240

Visual Inspection

Product should be clear and colorless, but it is acceptable for the reagent to contain a white flocculent material of a non-biological origin.

Bioburden

Visible bacterial colonies must be 100 or less CFU/mL


vincent A. Faddis
Quality Department Manager

10-Feb-2020
Date