

This document contains confidential, proprietary information of Ventana Medical Systems, Inc. It may not be copied or reproduced without prior written permission from Ventana Medical Systems, Inc.	Document Part No.	RFT-006-01	REV	Page
	RA Parent Doc No.	RI-006	A	1 of 1

TITLE: D134691 Rev. 2– Quality Control Specification

Product Name ULTRA Cell Conditioning Solution (CC1)
Ventana REF 950-224; 951-224
Roche Number 05424569001; 05424569214

Final Acceptance Testing pH acceptable range must be within 8.55 +/- 0.1 at 22°C
 Conductivity acceptable range must be within 0.22 mS/cm to 0.38 mS/cm.

Visual Inspection
 The reagent must be clear and free of turbidity.
 The reagent must be free of foreign material and substantially free of particulate matter.

Quality Department Manager

Date

< End of Document >