

This document contains confidential, proprietary information of Ventana Medical Systems, Inc. It may not be copied or reproduced without prior written permission from Ventana Medical Systems, Inc.	<b>Document Part No.</b>	<b>RFT-006-01</b>	<b>REV</b>	<b>Page</b>
	<b>RA Parent Doc No.</b>	<b>RI-006</b>	<b>A</b>	<b>1 of 1</b>

**TITLE: D134716 Rev. 3 – Quality Control Specification**

---

**Product Name** Cell Conditioning Solution (CC1)

**Ventana** REF 950-124; 951-124

**Roche Number** 05279801001; 05279801214

Final Acceptance Testing pH acceptable range must be within 8.55 +/- 0.1 at 22°C

Conductivity acceptable range must be within 0.22 mS/cm to 0.38 mS/cm.

Visual Inspection

The reagent must be clear and free of turbidity.

The reagent must be free of foreign material and substantially free of particulate matter.

---

**Quality Department Manager**

**Date**

**< End of Document >**