


TECHNICAL DATA SHEET

Product: Disposable High Pressure Syringe

REF# SMR101

Applicable injector: MEDRAD Vistron, Envision CT



No.	Component Name	Materials	No.	Component Name	Materials
①+②	Plunger assembly	ABS+Synthetic Rubber	⑥, ⑧	Luer lock	PC
③	Barrel (200ml)	PET	⑦	tube	PVC
④, ⑤, ⑨	cap	PP	⑩	Quick Fill Tube	PE

Sterilization method: EO Ethylene Oxide

Shelf life: 3 years after sterilization

Primary packaging: PS Blister tray sealed with Tyvek paper

Secondary packaging: Corrugated box

Pressure limit: 350psi

UNION MEDICAL SHENZHEN CO., LTD.

Add/ Room 603, Building 3, Fantasia MIC Plaza, Nanhai Ave., Nanshan District, 518062 Shenzhen, China

Tel/ +86 755 26696966 **Fax/** +86 755 26818966 **Web/** www.unionmed.net www.unionmedsyringe.com