

Technical Specification

Product name: CONCENTRATED SYSTEM LIQUID

Part Number: BO11524

1. Intended use: IVD Clinical Chemistry Analysers Consumable

2. Packaging Specification

Bottle: HDPE

Cap: PP

Packaging: Cartonbox

3. Physical and Chemical Specification

Storage: 15-30°C

Attribute	Test method	Specification
Physical State	NA	liquid
Colour	NA	colourless
Volume	NA	1 L
pH Test	NA	4.5
Absorbance	Prepare a 1/500 dilution in water (pipette 0.1 mL into a 50 mL volumetric flask) Measure the absorbance at 275 nm	0,300 - 0,500
Conductivity	NA	NA



BioSystems

Costa Brava, 30, 08030 BCN
NIF: A08678823

Laura López

Regulatory Affairs

Barcelona, 6th July 2022