

NANNING PASSION MEDICAL EQUIPMENT CO., LTD.

Quality Control Certificate

Product name:	MFG. Date:	Quantity:	Test Date:	
Electric chemotherapy chair PY-YD-310	1-Mar-2021	30 units	1-Mar-2021	
Items	Basic Requirements		Result	Judgement
1.Appearance	1.1 Weld should be smooth, cannot have burns, cracking leakage welding, surface is neat and beautiful, there is no front edge burr.		√	qualified
	1.2 Bedstead and panel installation is firm, cannot have obvious tilt distortion.		√	qualified
	1.3 The hinge, joint parts flexible and reliable operation.		√	qualified
	1.4 Two sidebars installation is firm, rotating flexible, reliable locking.		√	qualified
	1.5 The mat should be smooth and beautiful.		√	qualified
2. Technical Specification	2.1 Overall chair dimension: 1950mm L x 600 mm W x 560~820 mm H		√	qualified
	2.2 Back section adjusting: -10° ~75°		√	qualified
	2.3 Leg section adjusting: -80° ~10°		√	qualified
	2.4 Trendelenburg: -12°		√	qualified
	2.5 Frame: Q235 steel		√	qualified
	2.6 Cushion: PVC leather		√	qualified
	2.7 Weight of Chair: 96kgs		√	qualified
	2.8 Safe Maximum Load: 240kg		√	qualified
	2.9 Covering: ABS		√	qualified
	2.10 Motors: 3 Linak motors		√	qualified
	2.11 Power supply: AC220V 50Hz British Standard		√	qualified
3.Technical Indicator	3.1 Head pillow can be moved		√	qualified
	3.2 Linak motors for different position		√	qualified
	3.3 PVC for anti-fire		√	qualified
	3.4 Each accessory should be installed firmly, easy to remove (IV ROPE)		√	qualified
RESULT: The product according to the YY0571 medical equipment standard, determine the qualified. The result is accord with the contract require.				

note: "√" "Requirements" or "qualified", "x" Said "do not conform to the requirements" or "unqualified"

QA Inspector: Qi Quancheng

Check by: Fan Shu

