

96x140x16.5mm

金色代表燙銀

Product Illustration:



How to Use the Hearing Amplifier?

1. The hearing amplifier comes with 2 different ear bud sizes (large & small). Please select the size that best fits your ear.
2. Place the hearing amplifier into your ear with the amplifier's hook behind your ear.
3. Adjust the volume to a suitable level. If you hear a harsh sound, it means the volume setting is too high or the hearing amplifier is not placed in your ear properly. Adjust the position of the hearing amplifier or reduce the volume setting slightly to resolve this problem.
4. Take the hearing amplifier out of your ear when you finish using it and place it in the portable case. Ensure the hearing amplifier is in the OFF mode (Take out the battery) before placing it in the case.

ISO 13485



Made in Taiwan

Hearing Amplifier
BTE Type Battery Operated

Mimitakara® JAPAN

Mimitakara® JAPAN

Hearing Amplifier
BTE Type Battery Operated



Mimitakara® JAPAN
BTE Type Battery Operated

Hearing Amplifier
BTE Type Battery Operated



Nguyễn Hoàng Phương