



135x50x180mm



Item No.	000Hhearing P1
Max. Output Sound Level Gain	45 ± 5 dB
Continuous Use	Up to 10 hours per charge
Operating Temperature	+10°C ~ +40°C
Operating Humidity	30-85% RH
Storage Temperature	-10°C ~ +60°C
Storage Humidity	30-85% RH
Dimension	45 x 33 x 14 mm
Weight	21.4 g

Copyright © Digibone 2016 V.1.0



Mimitakara®

Digital Pocket Size Type

Rechargeable Hearing Aid

Patented audio technology
Charge by Micro USB
Bluetooth and app compatible
Control with smartphones



Digital 8 Channel
Suitable for both ear

ISO 13485 CE 2460
FDA Registration# 3005053820
Owner/Operator# 9071124

Mimitakara®
Rechargeable Hearing Aid
Digital Pocket Size Type

Mimitakara®
Rechargeable Hearing Aid
Digital Pocket Size Type